

APOBEC3C

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:20 - 1:50

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human APOBEC3C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

A3C; PBI; ARP5; ARDC2; ARDC4; APOBEC1L; bK150C2.3;

Catalog #: A5877

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27350

Isotype: IgG

Swiss Prot: Q9NRW3

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22.

Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control.

To place an order, please [Click HERE](#).