

KTN1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 156kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human KTN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

hj

Synonym:

CG1; KNT; MU-RMS-40.19;

Catalog #: A5879

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3895

Isotype: IgG

Swiss Prot: Q86UP2

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral membrane protein that is a member of the kinectin protein family. The encoded protein is primarily localized to the endoplasmic reticulum membrane. This protein binds kinesin and may be involved in intracellular organelle motility. This protein also binds translation elongation factor-delta and may be involved in the assembly of the elongation factor-1 complex. Alternate splicing results in multiple transcript variants of this gene.

To place an order, please [Click HERE](#).