

GLRX3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GLRX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

GRX3; GRX4; GLRX4; PICOT; TXNL2; TXNL3;

Catalog #: A5892

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10539

Isotype: IgG

Swiss Prot: O76003

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the glutaredoxin family. Glutaredoxins are oxidoreductase enzymes that reduce a variety of substrates using glutathione as a cofactor. The encoded protein binds to and modulates the function of protein kinase C theta. The encoded protein may also inhibit apoptosis and play a role in cellular growth, and the expression of this gene may be a marker for cancer. Pseudogenes of this gene are located on the short arm of chromosomes 6 and 9. Alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).