

CHD2

Reactivity: Human Mouse Rat

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:20 - 1:50

Calculated MW: 211kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CHD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fj

Synonym:

EEOC;

Catalog #: A5895

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1106

Isotype: IgG

Swiss Prot: O14647

Purity: Affinity purification

For research use only.

Background:

The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

To place an order, please [Click HERE](#).