

ABT1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ABT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

hABT1;

Catalog #: A5898

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29777

Isotype: IgG

Swiss Prot: Q9ULW3

Purity: Affinity purification

For research use only.

Background:

Basal transcription of genes by RNA polymerase II requires the interaction of TATA-binding protein (TBP) with the core region of class II promoters. Studies in mouse suggest that the protein encoded by this gene likely activates basal transcription from class II promoters by interaction with TBP and the class II promoter DNA.

To place an order, please [Click HERE](#).