

OXGR1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human OXGR1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GPR80; GPR99; P2Y15; P2RY15;

To place an order, please [Click HERE](#).

Catalog #: A5903

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27199

Isotype: IgG

Swiss Prot: Q96P68

Purity: Affinity purification

For research use only.