

NCL

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:76kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NCL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C23;

Catalog #:A5904

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4691

Isotype:IgG

Swiss Prot:P19338

Purity:Affinity purification

For research use only.

Background:

Nucleolin (NCL), a eukaryotic nucleolar phosphoprotein, is involved in the synthesis and maturation of ribosomes. It is located mainly in dense fibrillar regions of the nucleolus. Human NCL gene consists of 14 exons with 13 introns and spans approximately 11kb. The intron 11 of the NCL gene encodes a small nucleolar RNA, termed U20.

To place an order, please [Click HERE](#).