

RHOA

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:20 - 1:200 IF 1:50 - 1:200

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RHOA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ARHA; ARH12; RHO12; RHOH12;

Catalog #: A5909

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 387

Isotype: IgG

Swiss Prot: P61586

Purity: Affinity purification

For research use only.

Background:

RhoA is a small GTPase protein known to regulate the actin cytoskeleton in the formation of stress fibers. It acts upon two known effector proteins: Rho kinase (ROCK) and mdia (Diaphenous protein). RhoA is part of a larger family of related proteins known as the Ras superfamily.

To place an order, please [Click HERE](#).