

CSTF1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 48.3kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CSTF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

CstF-50; CstFp50;

Catalog #: A5915

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1477

Isotype: IgG

Swiss Prot: Q05048

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of three subunits which combine to form cleavage stimulation factor (CSTF). CSTF is involved in the polyadenylation and 3' end cleavage of pre-mRNAs. Similar to mammalian G protein beta subunits, this protein contains transducin-like repeats. Several transcript variants with different 5' UTR, but encoding the same protein, have been found for this gene.

To place an order, please [Click HERE](#).