

EIF1AX

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 16kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EIF1AX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1;

Catalog #: A5917

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1964

Isotype: IgG

Swiss Prot: P47813

Purity: Affinity purification

For research use only.

Background:

This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNAi and eIF3) to the 5' end of capped RNA.

To place an order, please [Click HERE](#).