

NHP2L1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:14kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NHP2L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

FA1; FA-1; NHPX; 15.5K; OTK27; SNU13; SSFA1; SPAG12; SNRNP15-5;

Catalog #:A5926

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4809

Isotype:IgG

Swiss Prot:P55769

Purity:Affinity purification

For research use only.

Background:

Originally named because of its sequence similarity to the *Saccharomyces cerevisiae* NHP2 (non-histone protein 2), this protein appears to be a highly conserved nuclear protein that is a component of the

To place an order, please [Click HERE](#).