

POLR2B

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 134kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human POLR2B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

eij

Synonym:

RPB2; POL2RB; hsRPB2; hRPB140

Catalog #: A5928

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5431

Isotype: IgG

Swiss Prot: P30876

Purity: Affinity purification

For research use only.

Background:

This gene encodes the second largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. This subunit, in combination with at least two other polymerase subunits, forms a structure within the polymerase that maintains contact in the active site of the enzyme between the DNA template and the newly synthesized RNA.

To place an order, please [Click HERE](#).