

PWP2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:102kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PWP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

UTP1; PWP2H; EHOC-17;

To place an order, please [Click HERE](#).

Catalog #:A5929

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5822

Isotype:IgG

Swiss Prot:Q15269

Purity:Affinity purification

For research use only.