

## RBM3

**Reactivity:** Human Mouse

**Tested applications:** WB

**Recommended Dilution:** WB 1:200 - 1:2000

**Calculated MW:** 17kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human RBM3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

RNPL; IS1-RNPL;

**Catalog #:** A5931

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5935

**Isotype:** IgG

**Swiss Prot:** P98179

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is a member of the glycine-rich RNA-binding protein family and encodes a protein with one RNA recognition motif (RRM) domain. Expression of this gene is induced by cold shock and low oxygen tension. A pseudogene exists on chromosome 1. Multiple alternatively spliced transcript variants that are predicted to encode different isoforms have been characterized although some of these variants fit nonsense-mediated decay (NMD) criteria.

*To place an order, please [Click HERE](#).*