

TRIM25

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 71kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human TRIM25

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1µl

Synonym:

EFP; Z147; RNF147; ZNF147;

Catalog #: A5940

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7706

Isotype: IgG

Swiss Prot: Q14258

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to the cytoplasm. The presence of potential DNA-binding and dimerization-transactivation domains suggests that this protein may act as a transcription factor, similar to several other members of the TRIM family. Expression of the gene is upregulated in response to estrogen, and it is thought to mediate estrogen actions in breast cancer as a primary response gene.

To place an order, please [Click HERE](#).