

## FXR1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 70kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human FXR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

bj

**Synonym:**

FXR1P;

**Catalog #:** A5942

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 8087

**Isotype:** IgG

**Swiss Prot:** P51114

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with polyribosomes, predominantly with the 60S ribosomal subunit. Three transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*