

RBM39

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human RBM39

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

HCC1; CAPER; RNPC2; FSAP59; CAPERalpha;

Catalog #: A5950

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9584

Isotype: IgG

Swiss Prot: Q14498

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the U2AF65 family of proteins. The encoded protein is found in the nucleus, where it co-localizes with core spliceosomal proteins. It has been shown to play a role in both steroid hormone receptor-mediated transcription and alternative splicing, and it is also a transcriptional coregulator of the viral oncoprotein v-Rel. Multiple transcript variants have been observed for this gene. A related pseudogene has been identified on chromosome X.

To place an order, please [Click HERE](#).