

RBM6

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 129kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human RBM6

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

3G2; g16; DEF3; DEF-3; HLC-11; NY-LU-12;

To place an order, please [Click HERE](#).

Catalog #: A5954

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 10180

Isotype: IgG

Swiss Prot: P78332

Purity: Affinity purification

For research use only.