

CPSF6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:3000

Calculated MW: 59kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CPSF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CFIM; CFIM68; HPBR11-4; HPBR11-7;

Catalog #: A5963

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11052

Isotype: IgG

Swiss Prot: Q16630

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. The cleavage factor complex is composed of four polypeptides. This gene encodes the 68kD subunit. It has a domain organization reminiscent of spliceosomal proteins.

To place an order, please [Click HERE](#).