

STRAP

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human STRAP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

MAWD; PT-WD; UNRIP;

To place an order, please [Click HERE](#).

Catalog #: A5964

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11171

Isotype: IgG

Swiss Prot: Q9Y3F4

Purity: Affinity purification

For research use only.