

CASC3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 76kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CASC3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

BTZ; MLN51;

Catalog #: A5967

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22794

Isotype: IgG

Swiss Prot: O15234

Purity: Affinity purification

For research use only.

Background:

The product of this gene is a core component of the exon junction complex (EJC), a protein complex that is deposited on spliced mRNAs at exon-exon junctions and functions in nonsense-mediated mRNA decay (NMD). The encoded protein binds RNA and interacts with two other EJC core components. It is predominantly located in the cytoplasm, but shuttles into the nucleus where it localizes to nuclear speckles.

To place an order, please [Click HERE](#).