

ERAL1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ERAL1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ERA; H-ERA; ERAL1A; HERA-A; HERA-B;

To place an order, please [Click HERE](#).

Catalog #: A5974

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26284

Isotype: IgG

Swiss Prot: O75616

Purity: Affinity purification

For research use only.