

HTATSF1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 86kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HTATSF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

TATSF1; TAT-SF1; dJ196E23.2;

Catalog #: A5977

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27336

Isotype: IgG

Swiss Prot: O43719

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene functions as a cofactor for the stimulation of transcriptional elongation by HIV-1 Tat, which binds to the HIV-1 promoter through Tat-TAR interaction. This protein may also serve as a dual-function factor to couple transcription and splicing and to facilitate their reciprocal activation. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).