

LARP7

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 67kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human LARP7

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ALAZS; PIP7S; HDCMA18P;

Catalog #: A5983

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51574

Isotype: IgG

Swiss Prot: Q4G0J3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is found in the 7SK snRNP (small nuclear ribonucleoprotein).

This snRNP complex inhibits a cyclin-dependent kinase, positive transcription elongation factor b, which is required for paused RNA polymerase II at a promoter to begin transcription elongation. A pseudogene of this gene is located on chromosome 3. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).