

RBM26

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 111kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human RBM26

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ARRS2; SE70-2; ZC3H17; PRO1777; C13orf10;

To place an order, please [Click HERE](#).

Catalog #: A5997

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64062

Isotype: IgG

Swiss Prot: Q5T8P6

Purity: Affinity purification

For research use only.