

GRWD1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 49kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human GRWD1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

CDW4; GRWD; RRB1; WDR28;

Catalog #: A6002

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83743

Isotype: IgG

Swiss Prot: Q9BQ67

Purity: Affinity purification

For research use only.

Background:

This gene encodes a glutamate-rich protein that contains five WD-repeat motifs. The encoded protein may play a critical role in ribosome biogenesis and may also play a role in histone methylation through interactions with CUL4-DDB1 ubiquitin E3 ligase.

To place an order, please [Click HERE](#).