

AQR

Reactivity:Human Mouse

Tested applications:WB ICC IP

Recommended Dilution:WB 1:500 - 1:1000 ICC 1:50 - 1:200 IP 1:20 - 1:50

Calculated MW:171kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human AQR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

IBP160; fSAP164;

To place an order, please [Click HERE](#).

Catalog #:A6011

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9716

Isotype:IgG

Swiss Prot:O60306

Purity:Affinity purification

For research use only.