

SRSF10

Reactivity: Human Mouse Hamster

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:1000 IP 1:20 - 1:50

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SRSF10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NSSR; TASR; SRp38; TASR1; TASR2; FUSIP1; FUSIP2; SFRS13; SRp40; SFRS13A; PPP1R149;

Catalog #: A6024

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10772

Isotype: IgG

Swiss Prot: O75494

Purity: Affinity purification

For research use only.

Background:

This gene product is a member of the serine-arginine (SR) family of proteins, which is involved in constitutive and regulated RNA splicing. Members of this family are characterized by N-terminal RNP1 and RNP2 motifs, which are required for binding to RNA, and multiple C-terminal SR/RS repeats, which are important in mediating association with other cellular proteins. This protein can influence splice site selection of adenovirus E1A pre-mRNA. It interacts with the oncoprotein TLS, and abrogates the influence of TLS on E1A pre-mRNA splicing. This gene has multiple pseudogenes. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms. In addition, transcript variants utilizing alternative polyA sites exist.

To place an order, please [Click HERE](#).