

GRSF1

Reactivity: Human

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:1000 IP 1:20 - 1:50

Calculated MW: 53kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human GRSF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Catalog #: A6027

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2926

Isotype: IgG

Swiss Prot: Q12849

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cellular protein that binds RNAs containing the G-rich element. The protein is localized in the cytoplasm, and has been shown to stimulate translation of viral mRNAs in vitro. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).