

HDLBP

Reactivity:Human Mouse

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:1000 IP 1:20 - 1:50

Calculated MW:141kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human HDLBP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

HBP; VGL; PRO2900;

Catalog #:A6028

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3069

Isotype:IgG

Swiss Prot:Q00341

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene binds high density lipoprotein (HDL) and may function to regulate excess cholesterol levels in cells. The encoded protein also binds RNA and can induce heterochromatin formation. Three transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).