

MAGOH

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:17kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant Protein of human MAGOH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAGOH1; MAGOHA;

Catalog #:A6035

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4116

Isotype:IgG

Swiss Prot:P61326

Purity:Affinity purification

For research use only.

Background:

Drosophila that have mutations in their mago nashi (grandchildless) gene produce progeny with defects in germplasm assembly and germline development. This gene encodes the mammalian mago nashi homolog. In mammals, mRNA expression is not limited to the germ plasm, but is expressed ubiquitously in adult tissues and can be induced by serum stimulation of quiescent fibroblasts.

To place an order, please [Click HERE](#).