

NOVA2

Reactivity: Human Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 49kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human NOVA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ANOVA; NOVA3;

To place an order, please [Click HERE](#).

Catalog #: A6039

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4858

Isotype: IgG

Swiss Prot: Q9UNW9

Purity: Affinity purification

For research use only.