

PIWIL2

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:3000

Calculated MW: 109kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human PIWIL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CT80; HILI; mili; PIWIL1L;

Catalog #: A6044

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55124

Isotype: IgG

Swiss Prot: Q8TC59

Purity: Affinity purification

For research use only.

Background:

PIWIL2 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells (Sasaki et al., 2003)

To place an order, please [Click HERE](#).