

PPP1R10

Reactivity: Human Mouse Hamster

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:1000 IP 1:20 - 1:50

Calculated MW: 99kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human PPP1R10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p99; FB19; R111; CAT53; PNUTS; PP1R10;

Catalog #: A6049

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5514

Isotype: IgG

Swiss Prot: Q96QC0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein phosphatase 1 binding protein. The encoded protein plays a role in many cellular processes including cell cycle progression, DNA repair and apoptosis by regulating the activity of protein phosphatase 1. This gene lies within the major histocompatibility complex class I region on chromosome 6, and alternatively spliced transcript variants have been observed for this gene.

To place an order, please [Click HERE](#).