

PRPF4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:100

Calculated MW: 58kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRPF4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PRP4; HPRP4; Prp4p; HPRP4P; SNRNP60;

Catalog #: A6052

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9128

Isotype: IgG

Swiss Prot: O43172

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is part of a heteromeric complex that binds U4, U5, and U6 small nuclear RNAs and is involved in pre-mRNA splicing. The encoded protein also is a mitotic checkpoint protein and a regulator of chemoresistance in human ovarian cancer. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).