

PTBP2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:4000

Calculated MW: 57kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PTBP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

PTB; nPTB; PTBLP; brPTB; nPTB5; nPTB6; nPTB7; nPTB8;

Catalog #: A6054

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 58155

Isotype: IgG

Swiss Prot: Q9UKA9

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene binds to the intronic cluster of RNA regulatory elements, downstream control sequence (DCS). It is implicated in controlling the assembly of other splicing-regulatory proteins. This protein is very similar to the polypyrimidine tract binding protein but it is expressed primarily in the brain.

To place an order, please [Click HERE](#).