

SERBP1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SERBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

CGI-55; CHD3IP; HABP4L; PAIRBP1; PAI-RBP1;

To place an order, please [Click HERE](#).

Catalog #: A6060

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26135

Isotype: IgG

Swiss Prot: Q8NC51

Purity: Affinity purification

For research use only.