

SNRNP40

Reactivity: Human Mouse Rat Hamster

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:1000 IP 1:20 - 1:50

Calculated MW: 39kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SNRNP40

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

40K; SPF38; WDR57; PRP8BP; HPRP8BP; PRPF8BP;

Catalog #: A6064

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9410

Isotype: IgG

Swiss Prot: Q96DI7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of the U5 small nuclear ribonucleoprotein (snRNP) particle. The U5 snRNP is part of the spliceosome, a multiprotein complex that catalyzes the removal of introns from pre-messenger RNAs.

To place an order, please [Click HERE](#).