

SNRNP70

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SNRNP70

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

RPU1; Snp1; U1AP; U170K; U1RNP; RNPU1Z; SNRP70; U1-70K;

To place an order, please [Click HERE](#).

Catalog #: A6065

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6625

Isotype: IgG

Swiss Prot: P08621

Purity: Affinity purification

For research use only.