

## SRRM1

---

**Reactivity:** Human Mouse Hamster

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:1000

**Calculated MW:** 102kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human SRRM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

160-KD; POP101; SRM160;

*To place an order, please [Click HERE](#).*

**Catalog #:** A6066

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 10250

**Isotype:** IgG

**Swiss Prot:** Q8IYB3

**Purity:** Affinity purification

For research use only.