

TIAL1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TIAL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

TCBP; TIAR;

Catalog #: A6075

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7073

Isotype: IgG

Swiss Prot: Q01085

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of a family of RNA-binding proteins, has three RNA recognition motifs (RRMs), and binds adenine and uridine-rich elements in mRNA and pre-mRNAs of a wide range of genes. It regulates various activities including translational control, splicing and apoptosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. The different isoforms have been shown to function differently with respect to post-transcriptional silencing.

To place an order, please [Click HERE](#).