

KHDRBS2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:100 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human KHDRBS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SLM1; SLM-1; bA535F17.1;

To place an order, please [Click HERE](#).

Catalog #: A6102

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 202559

Isotype: IgG

Swiss Prot: Q5VWX1

Purity: Affinity purification

For research use only.