

SCNN1G

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 74kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SCNN1G

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PHA1; BESC3; ENaCg; SCNEG; ENaCgamma;

Catalog #: A6126

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6340

Isotype: IgG

Swiss Prot: P51170

Purity: Affinity purification

For research use only.

Background:

Nonvoltage-gated, amiloride-sensitive, sodium channels control fluid and electrolyte transport across epithelia in many organs. These channels are heteromeric complexes consisting of 3 subunits: alpha, beta, and gamma. This gene encodes the gamma subunit, and mutations in this gene have been associated with Liddle syndrome.

To place an order, please [Click HERE](#).