

SMARCC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP CHIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:200

ChIP 1:20 - 1:100

Calculated MW: 122kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SMARCC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

Rsc8; SRG3; SWI3; BAF155; CRACC1;

Catalog #: A6128

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6599

Isotype: IgG

Swiss Prot: Q92922

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and contains a predicted leucine zipper motif typical of many transcription factors.

To place an order, please [Click HERE](#).