

PIGR

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 83kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PIGR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6130

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5284

Isotype: IgG

Swiss Prot: P01833

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the immunoglobulin superfamily. The encoded poly-Ig receptor binds polymeric immunoglobulin molecules at the basolateral surface of epithelial cells; the complex is then transported across the cell to be secreted at the apical surface. A significant association was found between immunoglobulin A nephropathy and several SNPs in this gene.

To place an order, please [Click HERE](#).