

AZU1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human AZU1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AZU; HBP; NAZC; AZAMP; CAP37; HUMAZUR;

Catalog #: A6132

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 566

Isotype: IgG

Swiss Prot: P20160

Purity: Affinity purification

For research use only.

Background:

Azurophil granules, specialized lysosomes of the neutrophil, contain at least 10 proteins implicated in the killing of microorganisms. The protein encoded by this gene is an azurophil granule antibiotic protein, with monocyte chemotactic and antimicrobial activity. It is also an important multifunctional inflammatory mediator. This encoded protein is a member of the serine protease gene family but it is not a serine proteinase, because the active site serine and histidine residues are replaced. The genes encoding this protein, neutrophil elastase 2, and proteinase 3 are in a cluster located at chromosome 19pter. All 3 genes are expressed coordinately and their protein products are packaged together into azurophil granules during neutrophil differentiation.

To place an order, please [Click HERE](#).