

TRPM2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 171kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRPM2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KNP3; EREG1; TRPC7; LTRPC2; NUDT9H; NUDT9L1;

Catalog #: A6137

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7226

Isotype: IgG

Swiss Prot: O94759

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a calcium-permeable cation channel that is regulated by free intracellular ADP-ribose. The encoded protein is activated by oxidative stress and confers susceptibility to cell death. Several alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known.

To place an order, please [Click HERE](#).