

IL31RA

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 82kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IL31RA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384;

Catalog #: A6141

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 133396

Isotype: IgG

Swiss Prot: Q8NI17

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the type I cytokine receptor family. This receptor, with homology to gp130, is expressed on monocytes, and is involved in IL-31 signaling via activation of STAT-3 and STAT-5. It functions either as a monomer, or as part of a receptor complex with oncostatin M receptor (OSMR). Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

To place an order, please [Click HERE](#).