

KMT2E

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:100

Calculated MW: 205kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human KMT2E

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fk

Synonym:

MLL5; NKp44L; HDCMC04P;

Catalog #: A6142

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55904

Isotype: IgG

Swiss Prot: Q8IZD2

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the myeloid/lymphoid or mixed-lineage leukemia (MLL) family and encodes a protein with an N-terminal PHD zinc finger and a central SET domain. Overexpression of the protein inhibits cell cycle progression. Alternate transcriptional splice variants have been characterized.

To place an order, please [Click HERE](#).