

SETDB1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 143kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SETDB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ESET; KG1T; KMT1E; TDRD21; H3-K9-HMTase4;

Catalog #: A6145

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9869

Isotype: IgG

Swiss Prot: Q15047

Purity: Affinity purification

For research use only.

Background:

This gene encodes a histone methyltransferase which regulates histone methylation, gene silencing, and transcriptional repression. This gene has been identified as a target for treatment in Huntington Disease, given that gene silencing and transcription dysfunction likely play a role in the disease pathogenesis. Alternatively spliced transcript variants of this gene have been described.

To place an order, please [Click HERE](#).